



Instructions for Use

Metanephrine ELISA


Enzyme immunoassay for the quantitative determination of
metanephrine in human urine

REF

40-371-25009



96



2-8°C

GenWay Biotech Inc.

6777 Nancy Ridge Drive • San Diego • CA 92121
Phone: 858-458-0866
Fax: 858-458-0833

techline@genwaybio.com
www.genwaybio.com

1. INTENDED USE

Enzyme immunoassay for the manual and automated quantitative determination of metanephrine in human urine.

2. SUMMARY AND EXPLANATION

The catecholamines adrenalin, noradrenalin and dopamine are synthesized in the adrenal medulla, the sympathetic nervous system and in the brain. They influence virtually all tissues and are involved together with other hormonal and neuronal systems in the regulation of a wide variety of physiological processes. As catecholamines and their metabolites metanephrine and normetanephrine are secreted in increasing amounts in a number of diseases.

In this context, diagnosis as well as the follow-up of tumor diseases of the nervous system are of special importance. This applies primarily to the pheochromocytoma, but also the neuroblastoma and the ganglioneuroma.

Malignant growth is described in 10% of pheochromocytomas. Furthermore, an increase of catecholamines and their metabolites metanephrine and normetanephrine can be observed in the carcinoid.

3. TEST PRINCIPLE

The assay procedure follows the basic principle of competitive ELISA whereby there is competition between a biotinylated and a non-biotinylated antigen for a fixed number of antibody binding sites. The amount of biotinylated antigen bound to the antibody is inversely proportional to the analyte concentration of the sample. When the system is in equilibrium, the free biotinylated antigen is removed by a washing step and the antibody bound biotinylated antigen is determined by use of anti-biotin alkaline phosphatase as marker and p-nitrophenyl phosphate as substrate. Quantification of unknowns is achieved by comparing the enzymatic activity of unknowns with a response curve prepared by using known standards.

4. WARNINGS AND PRECAUTIONS


1. For research use only. For professional use only.
2. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood. For further information (clinical background, test performance, automation protocols, alternative applications, literature, etc.) please refer to the GenWay Biotech Inc.-Homepage.
3. In case of severe damage of the kit package please contact GenWay Biotech Inc. in written form, latest one week after receiving the kit. Do not use damaged components in test runs, but keep safe for complaint related issues.
4. Obey lot number and expiry date. Do not mix reagents of different lots. Do not use expired reagents.
5. Follow good laboratory practice and safety guidelines. Wear lab coats, disposable latex gloves and protective glasses where necessary.
6. Reagents of this kit containing hazardous material may cause eye and skin irritations. See MATERIALS SUPPLIED and labels for details. Material Safety Data Sheets for this product are available on the GenWay Biotech Inc.-Homepage or upon request directly from GenWay Biotech Inc..
7. Chemicals and prepared or used reagents have to be treated as hazardous waste according to national biohazard and safety guidelines or regulations.
8. Avoid contact with Stop solution. It may cause skin irritations and burns.

5. STORAGE AND STABILITY

The kit is shipped at ambient temperature and should be stored at 2-8°C. Keep away from heat or direct sun light. The storage and stability of specimen and prepared reagents is stated in the corresponding chapters.

The microtiter strips are stable up to the expiry date of the kit in the broken, but tightly closed bag when stored at 2-8°C.

6. SPECIMEN COLLECTION AND STORAGE


	The in-vivo catecholamine and metanephrines release is influenced by several foods and drugs. Vitamin B, coffee and bananas, alpha-methyldopa, MAO and COMT inhibitors as well as medications related to hypertension should be discontinued for at least 72 h prior to specimen collection.
---	--

Urine

It is possible to use spontaneous as well as 24 h urine. The total volume of urine excreted during a 24 h period should be collected and mixed in a single bottle containing 10 - 15 mL of 6 N HCl as preservative. Determine total volume for calculation of results. **Mix and centrifuge samples before use in the assay.**

Storage:	≤ -20°C (Aliquots)	Keep away from heat or direct sun light. Avoid repeated freeze-thaw cycles.
Stability:	6 mon	

7. MATERIALS SUPPLIED

	The reagents provided with this kit are sufficient for up to 96 single determinations, or up to 48 duplicates, Metanephrine and Normetanephrine each.
---	---

Quantity	Symbol	Component
1 x 50 mL	ASSAYBUF CONC	Assay Buffer, Concentrate (10x) Contains: phosphate buffer, BSA, 1 % NaN ₃ .
1 x 2.25 mL	ACYLREAG	Acylation Reagent Ready to use. Contains: dimethylformamide.
1 x 10 mL	INDICATORBUF	Indicator Buffer Purple colored. Ready to use. Contains: Tris buffer, phenol red (color change at pH < 7.5).
2 x 50 x	HYDRTUB	Hydrolyzation Tubes Disposable polystyrene tubes (uncoated). Additional Hydrolyzation Tubes are available under REF KEZZ 661.
1 x 20 mL	HCL	HCl Ready to use. 0.1 M HCl.
1 x 12x8	MTP	Microtiter Plate Break apart strips. Coated with anti-rabbit IgG (goat, polyclonal).
1 x 7 x 0.35 mL	CAL A-G	Standard A-G 0; 26; 64; 160; 400; 1000; 2500 µg/L 0; 0.13; 0.33; 0.81; 2.03; 5.08; 12.7 µmol/L Ready to use. Contains: Metanephrine, 0.1 M HCl.
1 x 2 x 0.5 mL	CONTROL 1+2	Control 1+2 Ready to use. Exact concentrations see vial labels or QC Certificate.
3 x	BIOTIN LYO	Metanephrine Biotin, lyophilized Contains: Metanephrine Biotin, Tris buffer, < 0.1 % NaN ₃ (reconstituted).
1 x 7 mL	ANTISERUM	Metanephrine Antiserum Green colored. Ready to use. Contains: Antiserum (rabbit), phosphate buffer, 0.1 % NaN ₃ .
1 x 250 µL	ENZCONJ CONC	Enzyme Conjugate, Concentrate (100x) Contains: anti-Biotin antibodies, conjugated to alkaline phosphatase, Tris buffer, 0.01 % NaN ₃ .
9 x	PNPP SUBS	PNPP Substrate Tablets In one foil packet. Contains: p-nitrophenyl phosphate (PNPP).
1 x 27 mL	PNPP BUF	PNPP Substrate Buffer Ready to use. Contains: diethanolamine, water, 0.05 % NaN ₃ .
1 x 15 mL	PNPP STOP	PNPP Stop Solution Ready to use. Contains: 1 M NaOH, 0.25 M EDTA.
1 x 50 mL	WASHBUF CONC	Wash Buffer, Concentrate (10x) Contains: Tris buffer, HCl, Tween, 0.2 % NaN ₃ .
3 x	FOIL	Adhesive Foil

8. MATERIALS REQUIRED BUT NOT SUPPLIED

1. Micropipettes (Multipette Eppendorf or similar devices, < 3% CV). Volumes: 10; 50; 100; 1000 µL
2. Disposable glass test tubes (12 x 75 mm)
3. Orbital shaker (200-900 rpm) (e.g. EAS 2/4, SLT)
4. Vortex mixer
5. Water bath, 90°C
6. 8-Channel Micropipettor with reagent reservoirs
7. Wash bottle, automated or semi-automated microtiter plate washing system
8. Microtiter plate reader capable of reading absorbance at 405 nm (reference wavelength 600-650 nm)
9. Bidistilled or deionised water
10. Paper towels, pipette tips and timer

9. PROCEDURE NOTES

1. Any improper handling of samples or modification of the test procedure may influence the results. The indicated pipetting volumes, incubation times, temperatures and pretreatment steps have to be performed strictly according to the instructions. Use calibrated pipettes and devices only.
2. Once the test has been started, all steps should be completed without interruption. Make sure that required reagents, materials and devices are prepared ready at the appropriate time. Allow all reagents and specimens to reach room temperature (18-25 °C) and gently swirl each vial of liquid reagent and sample before use. Mix reagents without foaming.
3. Avoid contamination of reagents, pipettes and wells/tubes. Use new disposable plastic pipette tips for each component and specimen. Do not interchange caps. Always cap not used vials. Do not reuse wells/tubes or reagents.
4. It is advised to determine samples in duplicate to be able to identify potential pipetting errors.
5. Use a pipetting scheme to verify an appropriate plate layout.
6. Incubation time affects results. All wells should be handled in the same order and time sequences. It is recommended to use an 8-channel Micropipettor for pipetting of solutions in all wells.
7. Microplate washing is important. Improperly washed wells will give erroneous results. It is recommended to use a multichannel pipette or an automatic microplate washing system. Do not allow the wells to dry between incubations. Do not scratch coated wells during rinsing and aspiration. Rinse and fill all reagents with care. While rinsing, check that all wells are filled precisely with Wash Buffer, and that there are no residues in the wells.
8. Humidity affects the coated wells/tubes. Do not open the pouch until it reaches room temperature. Unused wells/tubes should be returned immediately to the resealed pouch including the desiccant.

10. PRE-TEST SETUP INSTRUCTIONS

- For manual and automatic version.



The contents of the kit for 96 determinations can be divided into 3 separate runs.

The volumes stated below are for one run with 4 strips (32 determinations).

If the customer wants to reduce the number of standards from 7 to 6 he can omit Standard G. The reportable range will then be reduced to 3000 µg/L.

If a larger number of strips is to be used, the volumes have to be changed accordingly.

10.1. Preparation of lyophilized or concentrated components




Do not mix up Metanephrine and Normetanephrine Enzyme Conjugate in case you use the Normetanephrine ELISA in parallel.

Dilute/dissolve	Component		Diluent	Relation	Remarks	Storage	Stability
15 mL	Assay Buffer	ad 150 mL	bidist. water	1:10		2-8°C	2 w
15 mL	Wash Buffer	ad 150 mL	bidist. water	1:10		2-8°C	4 w

1 vial	Metanephrine Biotin	with 2 mL	diluted Assay Buffer		Prepare freshly and use only once. Let stand for 15 min. Mix without foaming.	≤ -20°C (Aliquots)	2 mon
60 µL	Enzyme Conjugate	with 6 mL	diluted Assay Buffer	1:101	Prepare freshly and use only once.	18-25°C	5 h
3	PNPP Substrate Tablets	with 8 mL	PNPP Substrate Buffer		Prepare freshly and use only once.	18-25°C	5 h

10.2. Hydrolyzation of Urine Samples, Standards and Controls for total Metanephrine (in Hydrolyzation Tubes)

	<p>The hydrolyzation step is necessary for the determination of <u>total</u> normetanephrine and total metanephrine. No hydrolyzation is required when assaying <u>free</u> normetanephrine and metanephrine.</p> <p>Samples suspected to contain concentrations higher than the highest standard have to be diluted with 0.1 M HCl before hydrolyzation step.</p>
---	--

10.2.1. Sample preparation in the hydrolyzation tubes

1.	Pipette 10 µL of each Standard, Control and urine sample into the respective hydrolyzation tubes.
2.	Pipette 40 µL of 0.1 M HCl into each tube.
3.	Close tubes. Hydrolyze 1 h at 90°C (check temperature with thermometer). Allow to cool down to room temperature afterwards. Vortex.
4.	Pipette 100 µL of Indicator Buffer into each tube. Vortex.
5.	Pipette 20 µL of Acylation Reagent into each tube. Vortex each tube immediately after pipetting. Take care that addition of acylation reagent into the content of the tubes is complete.
6.	Close tubes. Incubate 15 min at Roomtemperature .
7.	Pipette 1 mL of diluted Assay Buffer into each tube. Vortex.

11. TEST PROCEDURE

In microtiter plate for manual and automated version

1.	Pipette 50 µL of each <u>acylated</u> Standard, acylated Control and <u>acylated</u> patient sample into the respective wells of the microtiter plate.
2.	Pipette 50 µL of Metanephrine Biotin into each well.
3.	Pipette 50 µL of Metanephrine Antiserum into each well.
4.	Cover plate with adhesive foil. Shake plate carefully. Incubate 1 h at RT on an orbital shaker (500 rpm) .
5.	Remove adhesive foil. Discard incubation solution. Wash plate with an automate 6 x with 250 µL of diluted Wash Buffer, (3 x manually) . Remove excess solution by tapping the inverted plate on a paper towel.
6.	Pipette 150 µL of freshly prepared Enzyme Conjugate into each well.
7.	Cover plate with new adhesive foil. Incubate 30 min at RT (18-25°C) on an orbital shaker (500 rpm).
8.	Remove adhesive foil. Discard incubation solution. Wash plate with an automate 6 x with 250 µL of diluted Wash Buffer, (3 x manually) . Remove excess solution by tapping the inverted plate on a paper towel.
9.	For adding of Substrate and Stop Solution use, if available, an 8-channel Micropipettor. Pipetting should be carried out in the same time intervals for Substrate and Stop Solution. Use positive displacement and avoid formation of air bubbles.
10.	Pipette 100 µL of freshly prepared PNPP Substrate Solution into each well.
11.	Incubate 40 min at RT (18-25°C) on an orbital shaker (500 rpm).
12.	Stop the substrate reaction by adding 100 µL of PNPP Stop Solution into each well. Briefly mix

	contents by gently shaking the plate.
13.	Measure optical density with a photometer at 405 nm (Reference-wavelength: 600-650 nm) within 60 min after pipetting of the Stop Solution.

12. QUALITY CONTROL

The test results are only valid if the test has been performed following the instructions. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable standards/laws. All kit controls must be found within the acceptable ranges as stated on the vial labels. If the criteria are not met, the run is not valid and should be repeated. Each laboratory should use known samples as further controls.

In case of any deviation the following technical issues should be proven: Expiration dates of (prepared) reagents, storage conditions, pipettes, devices, incubation conditions and washing methods.

It is recommended to participate at appropriate quality assessment trials.

13. CALCULATION OF RESULTS

The obtained OD of the standards (y-axis, linear) are plotted against their concentration (x-axis, logarithmic) either on semi-logarithmic graph paper or using an automated method. A good fit is provided with cubic spline, 4 Parameter Logisitcs or Logit-Log.

For the calculation of the standard curve, apply each signal of the standards (one obvious outlier of duplicates might be omitted and the more plausible single value might be used).

The concentration of the samples can be read from the standard curve.

In case of diluted samples the values have to be multiplied with the corresponding dilution factor.

Samples showing concentrations above the highest standard have to be diluted as described in PRE-TEST SETUP INSTRUCTIONS and reassayed.

Calculate the 24 h excretion for each urine sample: $\mu\text{g}/24\text{h} = \mu\text{g}/\text{L} \times \text{L}/24\text{h}$

Conversion:


1 ng/mL = 1 $\mu\text{g}/\text{L}$

Metanephrine ($\mu\text{g}/\text{L}$) $\times 5.07 \times 10^{-3} = \mu\text{mol}/\text{L}$

Typical Calibration Curve

(Example. Do not use for calculation!)

Standard	Metanephrine ($\mu\text{g}/\text{L}$)	Mean OD	OD/OD _{max} (%)
A	0.0	2.220	100.0
B	26	1.987	90.0
C	64	1.707	77
D	160	1.393	63
E	400	1.048	47
F	1000	0.702	32
G	2500	0.427	19



14. EXPECTED VALUES

Apparently healthy subjects show the following values:

Mean: 134 $\mu\text{g}/\text{d}$ Range: 25 – 312 $\mu\text{g}/\text{d}$ (95 % percentile)

It is recommended that each laboratory establishes its own range of normal values.

15. LIMITATIONS OF THE PROCEDURE

Specimen collection has a significant effect on the test results. See SPECIMEN COLLECTION AND STORAGE for details.

For cross-reactivities, see PERFORMANCE.

16. PERFORMANCE

Analytical Specificity (Cross Reactivity)	Substance	Cross Reactivity (%)		Cross-reactivity of other substances tested < 0.01 %	
	Adrenalin	0.327			
	Normetanephrine	0.151			
	Caffeic Acid	0.004			
Analytical Sensitivity (Limit of Detection)	8.6 µg/L	Mean signal (Zero-Standard) - 2SD			
Precision	Range (µg/L)	CV (%)			
	Intra-Assay	79.6 - 363	8.3 – 9.6		
	Inter-Assay	94.4 - 1650	8.1 – 9.3		
Linearity	Range (µg/L)	Serial dilution up to	Range (%)		
	94.6 - 6327	1:64	85 - 117		
Recovery	Mean (%)	Range (%)	% Recovery after spiking		
	102	85 - 114			
Method Comparison versus GCMS	GCMS = 1.04 x GenWay-Assay - 12.5			r = 0.88; n = 43	
Method comparison manual vs. automated version	Manual version = 0,97 Triturus + 0.380			r = 0.99; n = 56	