

Bead Name:	Anti-HSA IgY Coupled Beads
Lot Number:	LN 130619
<hr/>	
Formulation:	TBS (10 mM Tris pH 7.4, 150 mM NaCl, pH7.4) + 0.02% NaN ₃
Bead Concentration:	15 mg IgY/mL bead
Bead Volume :	0.5 mL
Storage Instructions:	4 °C, DO NOT FREEZE

Depletion Data:

One HSA Depletion Column can remove >99% of HSA from 30 - 40µl of human serum or plasma. The column can be re-used at least 20 times without a loss of binding efficiency.

Diluted serum or plasma before and after depletion were visualized on a SDS-PAGE gel.



M: MW marker, **IN:** Input, **F.T:** Flow Through depleted for HSA, **E:** Elution bound proteins from anti-HSA beads
4-20% SDS PAGE under reducing conditions.