



## Antibody Datasheet

**GenWay Biotech, Inc.**

6777 Nancy Ridge Drive

San Diego, CA 92121

Phone: 858.458.0866

Fax: 858.458.0833

Email: [techline@genwaybio.com](mailto:techline@genwaybio.com)

<http://www.genwaybio.com>

Rev: 042715

---

**Antibody Name:** Human Fibrinogen Antibody (Chicken IgY)  
**Description:** chicken polyclonal anti-Fibrinogen IgY antibody  
**Lot Number:** 140602

---

**Reactivity:** human Fibrinogen Protein  
**Immunogen:** human Fibrinogen Protein  
**Type:** polyclonal, IgY  
**Immunization host:** chicken  
**Purification:** Affinity purified

---

**Formulation:** 1x PBS Buffer, pH 7.4 + 0.02% NaN<sub>3</sub>  
**Concentration:** 2 mg/mL, by UV<sub>280</sub> nm.  
**Storage Instructions:** 4°C for short term (weeks), -20°C for long term. Aliquot to avoid frequent freeze/thaw.

---

**Applications:** Immunoassays

---

### QC Results:

Western Blot analysis using Human Fibrinogen Antibody (Chicken IgY) against Human Plasma sample under reducing conditions.

