



Antibody Datasheet

GenWay Biotech, Inc.

6777 Nancy Ridge Drive

San Diego, CA 92121

Phone: 858.458.0866

Fax: 858.458.0833

Email: techline@genwaybio.com

www.genwaybio.com

Bead Name: Anti-Transferrin IgY Coupled Beads

Lot Number: LN 150107

Formulation: TBS (10mM Tris-HCl, 150mM NaCl, pH 7.4)+ 0.02% NaN₃

Bead Concentration: 15 mg IgY/mL bead

Bead Volume : 0.5 mL

Storage Instructions: 4 °C, **DO NOT FREEZE**

Depletion Data:

One Transferrin Depletion Column can remove >99% of transferrin from 40 µl of human serum or plasma. One column can be re-used at least 20 times without a loss of binding efficiency.

Depletion Western Blot:

Diluted serum or plasma before and after depletion were visualized on the membrane from western blot. Transferrin bound to the anti-Transferrin IgY column was eluted and visualized by Western Blot.

